

# Crimpcap 11 mm



150016



150004



150007



150017



150019



150020

**Article**      **Description**      **11 mm crimpcaps**      **Packaging**

150016	Crimpcap, 11 mm, silver, 1.0 mm Natural Rubber/PTFE, red	Natural rubber/PTFE septa	1000
150004	Crimpcap, 11 mm, red, 1.0 mm, N.R./ butyl rubber/PTFE, red	Natural/butylrubber/PTFE septa	1000
150007	Crimpcap, 11 mm, gold, 1.0 mm, N.R./ butyl rubber/PTFE, red	Natural/butylrubber/PTFE septa	1000
150017	Crimpcap, 11 mm, silver, 1.0 mm, N.R./ butyl rubber/PTFE, red	Natural/butylrubber/PTFE septa	1000
150019	Crimpcap, 11 mm, green, 1.0 mm, N.R./ butyl rubber/PTFE, red	Natural/butylrubber/PTFE septa	1000
150020	Crimpcap, 11 mm, blue, 1.0 mm, N.R./ butyl rubber/PTFE, red	Natural/butylrubber/PTFE septa	1000



150022  
151130



150023  
151131



150062  
151133



151134



151132



970248

**Article**      **Description**      **11 mm crimpcaps**      **Packaging**

150022	Crimpcap, 11 mm, blank, 1.0 mm, Sil/PTFE, white/red	Sil/PTFE septa 1.0 mm	1000
150023	Crimpcap, 11 mm, green, 1.0 mm, Sil/PTFE, white/red	Sil/PTFE septa 1.0 mm	1000
150062	Crimpcap, 11 mm, red, 1.0 mm, Sil/PTFE, white/red	Sil/PTFE septa 1.0 mm	1000
151130	Crimpcap, 11 mm, blank, 1.2 mm Sil/PTFE, White/Red	Sil/PTFE septa 1.2 mm	1000
151131	Crimpcap, 11 mm, green, 1.2 mm Sil/PTFE, White/Red	Sil/PTFE septa 1.2 mm	1000
151132	Crimpcap, 11 mm, blue, 1.2 mm Sil/PTFE, White/Red	Sil/PTFE septa 1.2 mm	1000
151133	Crimpcap, 11 mm, red, 1.2 mm Sil/PTFE, White/Red	Sil/PTFE septa 1.2 mm	1000
151134	Crimpcap, 11 mm, gold, 1.2 mm Sil/PTFE, White/Red	Sil/PTFE septa 1.2 mm	1000
970248	Crimpcap, 11 mm, silver, 1.3 mm Sil/PTFE, White/Red, cross-slit 4 mm	Sil/PTFE septa 1.3 mm	1000



150012



151140

**Article**      **Description**      **11 mm crimpcaps**      **Packaging**

150012	Crimpcap, 11 mm, blank, 1.0 mm, PTFE/Sil/PTFE, Red/White/Red	PTFE/Sil/PTFE septa 1.0 mm	1000
151140	Crimpcap, 11 mm, blank, 0.1 mm PTFE, cap+rill-rand	PTFE septa 1.0 mm	1000