


Glass / Quartz microfiber filter

Glass fiber filters

Glass fiber filters are made of 100% borosilicate micro glass fibers. These filters offer an excellent level of small particles retention (up to 1µm) and have a large loading capacity.

They are suitable for microfiltration of air, gasses and liquids.

The filters are resistant to temperatures up to 500°C, have a good chemical resistance and are biological inert.



Filter selection:

GF/A	Suspended solids in waste water, air analysis
GF/B	Suspensions that are hard to filter (thicker filter)
GF/C	Suspended solids in waste water
GF/D	Pre-filter for membranes
GF/E	Collection of proteins
GF/H	Suspended solids in waste water
GF/6	Water and wastewater analysis (with binder, hydrophilic, up to 180°C)

Similar to:

MGA, GF50, GA-55, GMF1
MGB, GF51, GA-140, GMF2
MGC, GF52, GC-50, GMF3
MGD, GF53, GD-120, GMF4
MGF, GF55, GF-75, GMF5
934-AH, GF30, GS-25, GMF6

GF/A			GF/B		
Article	Description	Qty.	Article	Description	Qty.
890025	Filter, glass fiber, GF/A, 1.6 µm 25 mm	100	891025	Filter, glass fiber, GF/B, 1.0 µm 25 mm	50
890048	Filter, glass fiber, GF/A, 1.6 µm 47 mm	100	891047	Filter, glass fiber, GF/B, 1.0 µm 47 mm	50
890054	Filter, glass fiber, GF/A, 1.6 µm 55 mm	100	891055	Filter, glass fiber, GF/B, 1.0 µm 55 mm	50
890070	Filter, glass fiber, GF/A, 1.6 µm 70 mm	100	891070	Filter, glass fiber, GF/B, 1.0 µm 70 mm	50
890090	Filter, glass fiber, GF/A, 1.6 µm 90 mm	100	891090	Filter, glass fiber, GF/B, 1.0 µm 90 mm	50
900876	Filter, glass fiber, GF/A, 1.6 µm 110 mm	100	891110	Filter, glass fiber, GF/B, 1.0 µm 110 mm	25
890126	Filter, glass fiber, GF/A, 1.6 µm 125 mm	100	891125	Filter, glass fiber, GF/B, 1.0 µm 125 mm	25
890151	Filter, glass fiber, GF/A, 1.6 µm 150 mm	100	891150	Filter, glass fiber, GF/B, 1.0 µm 125 mm	25
GF/C			GF/D		
Article	Description	Qty.	Article	Description	Qty.
892025	Filter, glass fiber, GF/C, 1.2 µm 25 mm	100	897025	Filter, glass fiber, GF/D, 2.7 µm 25 mm	50
892047	Filter, glass fiber, GF/C, 1.2 µm 47 mm	100	897047	Filter, glass fiber, GF/D, 2.7 µm 47 mm	50
892055	Filter, glass fiber, GF/C, 1.2 µm 55 mm	100	897055	Filter, glass fiber, GF/D, 2.7 µm 55 mm	50
892070	Filter, glass fiber, GF/C, 1.2 µm 70 mm	100	897070	Filter, glass fiber, GF/D, 2.7 µm 70 mm	50
892090	Filter, glass fiber, GF/C, 1.2 µm 90 mm	100	897090	Filter, glass fiber, GF/D, 2.7 µm 90 mm	50
892110	Filter, glass fiber, GF/C, 1.2 µm 110 mm	100	897110	Filter, glass fiber, GF/D, 2.7 µm 110 mm	25
892125	Filter, glass fiber, GF/C, 1.2 µm 125 mm	100	897125	Filter, glass fiber, GF/D, 2.7 µm 125 mm	25
892150	Filter, glass fiber, GF/C, 1.2 µm 150 mm	100	897150	Filter, glass fiber, GF/D, 2.7 µm 150 mm	25

Glass / Quartz microfiber filter

GF/F

Article	Description	Qty.
894025	Filter, glass fiber, GF/F, 0.7 µm 25 mm	100
894047	Filter, glass fiber, GF/F, 0.7 µm 47 mm	100
894055	Filter, glass fiber, GF/F, 0.7 µm 55 mm	100
894070	Filter, glass fiber, GF/F, 0.7 µm 70 mm	100
894090	Filter, glass fiber, GF/F, 0.7 µm 90 mm	50
894110	Filter, glass fiber, GF/F, 0.7 µm 110 mm	50
894125	Filter, glass fiber, GF/F, 0.7 µm 125 mm	50
894150	Filter, glass fiber, GF/F, 0.7 µm 150 mm	50

GF/H

Article	Description	Qty.
895025	Filter, glass fiber, GF/H, 1.5 µm 25 mm	100
895047	Filter, glass fiber, GF/H, 1.5 µm 47 mm	100
895055	Filter, glass fiber, GF/H, 1.5 µm 55 mm	100
895070	Filter, glass fiber, GF/H, 1.5 µm 70 mm	100
895090	Filter, glass fiber, GF/H, 1.5 µm 90 mm	100
895110	Filter, glass fiber, GF/H, 1.5 µm 110 mm	100
895125	Filter, glass fiber, GF/H, 1.5 µm 125 mm	100
895150	Filter, glass fiber, GF/H, 1.5 µm 150 mm	100

GF/F

Article	Description	Qty.
893025	Filter, glass fiber, GF/6, 25 mm, with binder	50
893055	Filter, glass fiber, GF/6, 55 mm, with binder	50
893090	Filter, glass fiber, GF/6, 90 mm, with binder	50


Quartz microfiber filters

Quartz Microfiber filters are made of 100% quartz. These filters offer the same technical specifications as glass fiber filters, except for the higher temperature resistance (up to 900°C) and they contain a low amount of trace elements.

Applications: Air pollution analysis, emission control, exhaust gas control.

Quartz

Article	Description	Qty.
896025	Filter, quartz microfiber, 1.5 µm 25 mm	50
896047	Filter, quartz microfiber, 1.5 µm 47 mm	50
896055	Filter, quartz microfiber, 1.5 µm 55 mm	50
896070	Filter, quartz microfiber, 1.5 µm 70 mm	50
896090	Filter, quartz microfiber, 1.5 µm 90 mm	50
896110	Filter, quartz microfiber, 1.5 µm 110 mm	50
896125	Filter, quartz microfiber, 1.5 µm 125 mm	50
896150	Filter, quartz microfiber, 1.5 µm 150 mm	50



Note: If the filter you need is not listed, please feel free to contact us.